LayerVoyager

Contents

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LayerVoyager

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- Interface Description
- File 、 Tools 、 Settings 、 Help

Basic Operations


- Launch Program and Load Image
- Select Measurement or Annotation Tool
- Perform Measurement and View Results
- Save Image

Notes

- System Requirements
- Important Notes

Overview

Interface Description

File

1. Open Image
 - Supports opening image files in *.pkg format.
 - A single *.pkg file may contain multiple images.
2. Save Image
 - Save the currently displayed image.

Tools

1. Tools

| TOOL | DESCRIPTION |
|---------------------------|--|
| SELECT | Select objects drawn on the image (cursor changes to a hand icon) for adjustment. |
| LINE | Perform measurements using a straight line. Endpoints can be adjusted when the cursor becomes a hand icon. |
| RECTANGLE | Perform measurements using a rectangle. |
| CIRCLE | Perform measurements using a circle. |
| THREE-POINT CIRCLE | Measure a circle by selecting any two edge points first. Click and hold the third point to generate the circle, then drag to adjust position before releasing the mouse. |
| DIAMETER | Measure the diameter length of a circle. |
| POLYGON | Measure using a polygon. Double-click to finish drawing. |
| LABEL | Display text annotations on the image. |
| CALIBRATION | Draw a reference line according to a scale bar and set the actual length and unit. |

2. Shape / Area
 - Fill Shape :
 - Fill the interior area of rectangles, circles, three-point circles, and polygons.
 - Calculate Area :
 - Rectangle / Circle / Three-Point Circle / Polygon → Calculate enclosed area
 - Line / Diameter → Calculate length
 - Layer Synchronization :


- The software allows switching between different images.
- When enabled, drawn objects appear on all images at the same position.
- When disabled, objects appear only on the current image.

3. Color (RGB)

- Palette — Change measurement object line color
- Text Color — Change annotation text color
- Background Color — Change text background color

4. Text / Drawing

- Font Size — Text display size
- Line Width — Thickness of measurement lines
- Text Position Options


- Label Text — Used only with the Label Tool.

Settings

1. Image Rotation

- Rotates all loaded images, supporting only 90, 180, and 270 degree rotations.
- Measured items can also be moved, deleted, or left unchanged.
- Be advised that performing a rotation will reset the undo buffer.

2. Language

- Switch application language.
 - Currently supported:
 1. Traditional Chinese
 2. English
 3. Simplified Chinese

Help

1. Open the shortcut key reference document.

Basic Operations

1. Launch Program and Load Image



If one of the two warning windows appears during startup, it indicates that the software license has expired or is about to expire.

- Left Warning Window:
 - A warning message appears, but the software can still be used after confirmation.
- Right Warning Window:
 - The program will automatically close.
 - Please follow the instructions to download the latest version from the MA-tek official website.

After opening an image:



- Use the dropdown menu to switch images.
- Hold the Tab key and scroll the mouse wheel to switch images quickly.

2. Select Measurement or Annotation Tool

1. Select Calibration first.
2. Locate the scale bar in the image.
3. Draw a reference line with the same length.
4. Enter the actual value and confirm.

- After calibration, measurements can be performed directly on objects of interest.

3. Perform Measurement and View Results

- Measurement results are displayed in real time on the image.
- Measured objects can be:
 - dragged
 - modified
 - resized
 - recolored
 - adjusted in line thickness

⚠ Note:


Polygons can only be moved after creation (no vertex editing).

Polygon Completion

- Double-click to finish drawing (the final click also counts as a point).

Editing Measurement Parameters

- Select the object first.
- A dashed black selection box will appear around it.

4. Save Image

- The saved file corresponds to the currently displayed image.
- One image is saved at a time.
- Saved images may have reduced resolution.

Notes

✧ System Requirements

- 🔧 Operating System: Windows 10 / Windows 11
- 🔧 Recommended Hardware: 8 GB RAM or above
- 🔧 Display Recommendation: Full HD resolution or higher

✧ Important Notes

- Large images may impact system performance.
- If image display issues occur:
 - Click Clear Image, or
 - Delete the temp folder inside the loaded image directory.

The temp folder stores all measurement information.

⚠ Saving is not supported when the file path contains Chinese characters.